

 **EASY**<sup>®</sup> *Ambulatory* **2**



# Expand your diagnostic reach to hospital beds and patients' homes

Run a complete Type 2 sleep study anywhere with Easy® Ambulatory 2 PSG.

Our Ambulatory solution utilizes the same Easy III software as our clinical PSG system with flexible workspace configurations, montages and customizable reports.

Easy Ambulatory 2 PSG is composed of a recorder, battery pack, the 32-channel color-coded amplifier, respiratory belts, accessories and modules for nasal pressure and snoring, oximetry, body position and limb movement.

Add Q-Video Mobile™ 2 with automatic infrared sensors to capture video of your patient's unattended sleep study.



32-channel amplifier is color-coded and comes with a wearable pouch



Recorder and battery lets you record up to 48 hours



Q-Video Mobile 2 has infrared lights for low-light recording



Cadwell Industries, Inc.  
909 N. Kellogg St. · Kennewick, WA 99336  
(800) 245-3001 · +1 (509) 735-6481 ph · +1 (509) 783-6503 fax  
[www.cadwell.com](http://www.cadwell.com) · [info@cadwell.com](mailto:info@cadwell.com)

